

AL-ZX1 Auto Lift for Spinergy ZX-1

Harmar's unique AL-ZX1 Lift is specifically designed to easily transport the Spinergy ZX-1, the innovative power add-on device that turns a manual wheelchair into an easy-glide power chair.

The custom platform and hold-down arm automatically secures the ZX-1 safely for transport on the rear of a vehicle featuring a Class I or greater hitch. When empty, the platform can be raised up and out of the way, making this lift the perfect ZX-1 accessory.

Features & Benefits

- Automatic securement
- Lightweight, unobtrusive platform
- One-button convenience
- Easy load-on design

Specifications	AL-ZX1
Accommodates	Spinergy ZX-1 power add-on
Lifting Capacity	100 lbs. (45.4 kg)
Installed Weight	69 lbs.
Securement	Automatic hold-down arm
Hitch	Class I or greater
Power switch	Keyed
Controls	One-touch design
Options	Swing-Away
Warranty	3-year transferable

