

Date \_\_\_\_\_  
 Sales Representative \_\_\_\_\_  
 Dealer \_\_\_\_\_  
 Address \_\_\_\_\_  
 City \_\_\_\_\_ State \_\_\_\_\_ Zip \_\_\_\_\_  
 Phone \_\_\_\_\_  
 Fax \_\_\_\_\_  
 Email \_\_\_\_\_  
 Opportunity Reference / Project Name \_\_\_\_\_

Quote    Order  
 Account Number \_\_\_\_\_  
 Purchase Order Number \_\_\_\_\_  
 Shipping:  Pick-up    Fed Ex    Palletized (premium)  
 Ship To Address \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 Contact \_\_\_\_\_  
 Phone \_\_\_\_\_  
 Projected Close Date \_\_\_\_\_

### CONFIGURATION

QTY \_\_\_\_\_ SL600 Std. 15'6" Rail  
 QTY \_\_\_\_\_ SL600FR-LH Folding Rail, Left Hand, 15'6" Rail  
 QTY \_\_\_\_\_ SL600FR-RH Folding Rail, Right Hand, 15'6" Rail

### TRACK, PARTS & OPTIONS

#### Track


QTY \_\_\_\_\_ SL600T3 – 3' Mid Section Assembly w/brackets  
 QTY \_\_\_\_\_ SL600T5 – 5' Mid Section Assembly w/brackets  
 QTY \_\_\_\_\_ SL600T8 – 8' Mid Section Assembly w/brackets  
 [ ] SL600C- Custom Cut \_\_\_\_\_ Ft in length (including overrun)  
 single unit orders only.

Must include additional mid-section(s) above if total above 15'6".

#### Parts & Options

QTY \_\_\_\_\_ SL305 – Key Lock on Lift (chassis)  
 QTY \_\_\_\_\_ SL600SB1 – Additional Short Brackets  
 QTY \_\_\_\_\_ SL600TB1 – Additional Tall Brackets  
 QTY \_\_\_\_\_ SL600CC3 – 3' Charger Cord Extension  
 QTY \_\_\_\_\_ SL600CC8 – 8' Charger Cord Extension

### HOW TO MEASURE



NOTE: If A to H is < 13", must order tall bracket.

#### Measurements

\_\_\_\_\_ **A to B** - Inches from top landing to bottom floor

\_\_\_\_\_ **A to H** - Inches from top landing to nearest obstruction - for example, a door

**Rise** \_\_\_\_\_ Inches   **Tread** \_\_\_\_\_ Inches   **Width** \_\_\_\_\_ Inches

7" Overrun recommended if top landing is clear  
 (This is added to the A to B)

#### Notes

\_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_